

MAXIZAINO

Professional EMS backpack



Capacious professional backpack designed to accommodate all the necessary equipment used in ALS procedures.

Included Equipment		Specific Characteristics	
BG02007A	PRO B-LIFE - PVC AD. RES. W/OVERP.VALVE R-VALVE	Dual transport possibility thanks to the padded shoulder straps and carrying handle	
BS04011D	BURN S (BA) 10 X 10 CM - BURN STOP DRESSING	Detachable and adjustable shoulder straps with adjustable cross band to provide more stability during transport	
BS04021D	BURN S (BA) 20 X 20 CM - BURN STOP DRESSING	High visibility thanks to the reflective bands	
CB03030A	PHARMA TANK 1 - SEMI-RIGID PHIAL HOLDERS	2 capacious front pockets with zip closure	
DG00150A	YANKAUER PROBE	Detachable semirigid inner compartment equipped with zip closure and pocket with strap closure	
DG00212C	CATHETER CONTROLLED SUCTION CH/FG 12	In the main body of the backpack are present:	
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS	- two side pockets with zip closure	
DG01411C	DGX 411 - AN.SPHIGMO. 2TUBES W/AD. CUFF ORAN.	- one pocket with zip closure on the inner side of the front flap and a pocket with strap closure	
DG02110B	ECONOMICAL PHONENDOSCOPE	- adjustable straps to fix a cylinder up to 2l	
DG08062D	DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES	- strap bands for attachment of internal pockets	
DG80005A	IME INTUBATION KIT COMPLETE		
FG00605A	FORCEPS MAGILL ADULT		
FG00607A	MAGILL PAEDIATRIC - TUBE-GUIDING FORCEPS		
FG03308B	GMX 308 - SURGICAL FORCEPS		
FG08810A	MOUTH WEDGE EXCELL		
GM00363G	KLEMMER FORCEPS		
GM02215K	REEL OF PLASTER 1,25 cm x 5 m		
GM02225K	ADHESIVE BANDAGE ROLL FABRIC 2,5CMX5M		
GM02250K	CLOTH PLASTER 5CM X 5M		
GM02405K	GAUZE BENDAGE 5M x5CM		
GM02407K	GAUZE BANDAGE 5M X 7CM		
GM02410K	GAUZE BANDAGE 5M X 10CM		
GM03050K	STERILE DISPOSABLE SYRINGE 5ML		
GM03100K	DISPOSABLE STERILE SYRINGE 10 ML		
GM04020K	PLASTERS - 40 PIECES		
GM04250K	GMX 003 – DISINFECTANT 250 ml		
GM04300K	HYDROGEN PEROXIDE 10 VOL. 250 ml		
IF00100C	T.I.V. SPENCER AUTOMATIC TOURNIQUETTE		
IF00421A	COLD PACK SPENCER ICE TNT SINGLE USE BOX 25 Pcs		
RM10310A	GUEDEL AIRWAY KIT 6 SIZES		
RM10880A	CUFFED MASK PVC SZ.0 AIR CUSHION		
RM10893A	CUFFED MASK SIZE. 4 PVC CO. COD WITHOUT VALVE		
RM30300A	COLLIN - STAINLESS STEEL TONGUE FORCEPS		
SC75011E	TRANSPORT BAG FOR JET COMPACT		
SC75100D	JET COMPACT 300 D		
ST00630A	ISOTHERMIC SHEET ADULT SIZES 160X220 CM.		
TA09020A	FLEX MASK SET WITH ORANGE PVC CASE		
		Technical data	
		Dimensions	500x420x220 mm
		Detachable compartment dimensions	395x140x200 mm
		Weight of the empty backpack	3,19 kg
		Suction device technical data	
		Suction flow (free flow)	32l/min
		Vacuum range	From 50 to 550 mmHg
		IP classification	IP12
		Weight	1,6 kg
		Medical devices CE marked according to dir. 93/42/CE	
		MAXIZAINO COMPLETE ALS-PROFESSIONAL EMS BACKPACK	
		CB04023A	

Rev.0 (13/05/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it