


BOA BAG

Buoyant soft bag



Bag to store first aid devices particularly suitable for use in environments including water rescue.

Specific features

- Double possibility of transport thanks to the carrying handles and shoulder strap
- PVC fabric limits infiltration of liquids and allow easy cleaning
- Opening from the top with zip closure
- High visibility thanks to the reflective inserts
- The padded structure and the dimensions increase buoyancy
- Semirigid internal detachable compartment PVC coated outside, internally padded with separators. The position of the separators is adjustable thanks to the strap fixing
- Elastic bands on the top flap to store phials or accessories

Technical data

Dimensions	460x310x300 ± 10 mm
Dimensions detachable compartment	380x200x210 ± 10 mm
Upper opening dimensions	320x230 mm
Weight	1,95 ± 0,03 kg
Materials	PVC, Nylon

BOA BAG - BUOYANT SOFT BAG
CB04026A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it