

# MAXIZAINO

## Professional EMS backpack



Capacious professional backpack designed to accommodate all the necessary equipment used in ALS procedures.

### Specific Characteristics

- Dual transport possibility thanks to the padded shoulder straps and carrying handle
- Detachable and adjustable shoulder straps with adjustable cross band to provide more stability during transport
- High visibility thanks to the reflective bands
- 2 capacious front pockets with zip closure
- Detachable semirigid inner compartment equipped with zip closure and pocket with strap closure
- In the main body of the backpack are present:
  - two side pockets with zip closure
  - one pocket with zip closure on the inner side of the front flap and a pocket with strap closure
  - adjustable straps to fix a cylinder up to 2l
  - strap bands for attachment of internal pockets

### Technical data

Dimensions	500x420x220 mm
Detachable compartment dimensions	395x140x200 mm
Weight	3,19 kg

MAXIZAINO EMPTY - PROFESSIONAL EMS BACKPACK  
**CB04040A**

Rev.1 (26/05/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)