

STARTER KIT

Ambulance kit



A bag equipped with handles and shoulder strap , containing some of the essential accessories for first aid in an ambulance.

The devices contained inside are compliant with the reference European Directives, when applicable.

TECHNICAL DATA

External dimensions	450x330x340 mm
Bag material	Nylon
Colour	Blue
Equipment	
Clip Matrix	Eva and alluminium clipboard
Glovex	Latex gloves Size L 100pcs package
Pro Spike	Disposable container for used needles
GMX401	Disinfectant towel
ICE PACK	Instant Ice pack
OXYMASK 150	High concentration mask for oxygen therapy with reservoir, latex Free
FLEX MASK	Mouth to mouth mask with unidirectional disposable valve, latex Free
Roll Up CB04054K	Fleece blanket STARTER KIT AMBULANCE KIT FOR FIRST USE



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>