


BLITZ

The rescuer's backpack



Backpack characterized by reduced thickness and ideal to accommodate devices used by rescuers in first aid.

Specific Characteristics

- Dual transport possibility thanks to the padded shoulder straps and carrying handle
- Adjustable shoulder straps
- One capacious front pocket with zip closure
- In the main body of the backpack are present:
 - a mesh pocket on the frontal flap
 - 3 adjustable strap bands and elastic bands to store accessories
- The backpack can be totally opened, to allow an easy access to the content

Technical data

Dimensions	480x300x120 ± 20 mm
Volume	12 ± 1L
Material	Nylon, TPU
Weight empty	0,8 ± 0,1 kg

BLITZ EMPTY - THE RESCUER'S BACKPACK
CB04081K

Rev.1 (14/11/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it