

## R-AID Trolley

### Multipurpose trolley



Emergency trolley extremely versatile and robust. The system of 4 coaxial wheels ensures high load resistance while dragging on the ground.

#### Specific Characteristics

Excellent portability thanks to the trolley system, 2 handles and shoulder strap

1 main compartment with straps for fixing the tank (max 2l), and 2 transparent pockets with strap closure.

2 Side pockets with zip closure

1 Front transparent pocket with strap closure

1 Front pocket with zip closure

Front detachable backpack with padded shoulder straps

Reflective inserts for greater visibility

Reinforced bottom

#### Standard equipment

4 R-Series internal bags with strap fixing

#### Technical data

Dimensions	400x300x570 mm
Weight	3,8 kg
Volume	40 l
Materials	Nylon, PVC, PE, AL
Maximum Load	30 kg

R-AID TROLLEY – MULTIPURPOSE TROLLEY  
RED/BLACK W/BAGS  
**CB04130B**

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)