

R-AID Trolley

Multipurpose trolley



Emergency trolley extremely versatile and robust. The system of 4 coaxial wheels ensures high load resistance while dragging on the ground.

Specific Characteristics

Excellent portability thanks to the trolley system, 2 handles and shoulder strap

1 main compartment with straps for fixing the tank (max 2l), and 2 transparent pockets with strap closure.

2 Side pockets with zip closure

1 Front transparent pocket with strap closure

1 Front pocket with zip closure

Front detachable backpack with padded shoulder straps

Reflective inserts for greater visibility

Reinforced bottom

Technical data

Dimensions	400x300x570 mm
Weight	2,7 kg
Volume	40 l
Materials	Nylon, PVC, PE, AL
Maximum Load	30 kg

R-AID TROLLEY – MULTIPURPOSE TROLLEY

RED/BLACK EMPTY

CB04132B

Accessories

R-Series bags

CB90003A	R1 bag – Suitable for intubation
CB90005A	R2 bag – Suitable for ventilation
CB90002A	R3 bag – Suitable for medicines
CB90004A	R4 bag – Suitable for infusion

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it