

# R-AID Trolley PRO

## Multipurpose backpack with trolley



Emergency trolley extremely versatile and robust, usable also as a backpack. 4 coaxial wheels ensure high load resistance while dragging the trolley on the ground.

### Specific Characteristics

Excellent portability thanks to the trolley system, 2 handles and shoulder strap

1 main padded compartment, equipped with straps to fix a bottle (max 2l) and removable internal dividers that allow to split the compartment up to 9 parts

1 transparent pocket with zip closure on the internal side of the closure flap of the main compartment

2 Transparent pockets with zip closure of which:

- 1 equipped with internal pocket with zip closure
- 1 equipped with *transparent detachable bag* with zip closure, fixed by means of straps

1 Front transparent pocket with strap closure

1 Front pocket with zip closure

1 top pocket with zip closure and carrying handle

Padded shoulder straps and lumbar belt with quick release buckle improve load distribution

Reflective inserts for increased visibility

Reinforced bottom

### Equipment

5 M-Series internal bags with strap fixing

### Technical data

Dimensions	620x350x410 ± 10 mm
PVC bag dimensions	360x160 mm
Weight <sup>1</sup> (without internal bags)	3 ± 0,2 kg
Volume <sup>1</sup>	48 ± 0,5 l
Shoulder strap length	From 70 to 130 cm
Trolley handle extension	43 cm
Materials	Nylon, PVC, AI
Loading capacity	30 kg

<sup>1</sup> Main compartment volume 32 l

R-AID TROLLEY PRO - RED/BLACK CORDURA TROLLEY W/BAGS  
**CB04181B**

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)