

R-AID PRO

Multipurpose backpack



Multipurpose emergency backpack extremely versatile and robust. It is characterized by a structure integrated with the carrying points giving added strength and for a better support of the weight.

Specific features

Cordura fabric
Carrying point fixed to the supporting structure , with the internal part made of metal alloy and rigid polyethylene
TPE reinforcement on the bottom
2 external pockets with zip closure
1 removable top pocket with zip closure and equipped with elastic rings and pockets
1 removable front backpack attached by zip equipped with various straps and elastic rings, two shoulder straps and transparent front pocket
The padded main compartment provides maximum protection to the content and allows to store a set of bags
Strap bands to fix pockets or bags
Belts to fix an oxygen tank of 2lt
Ergonomic handles with loop for hanging kit and adjustable upper carrying point
Padded shoulder straps adjustable also in height in order to fit each operator
The padded lumbar support has a protection for the shoulder straps in order to protect them from dirt or damage
Anatomic padded back to provide maximum comfort
Reflective inserts increase visibility

Equipment

5 M series bags

Technical data

Dimensions	600x420x300 ± 30 mm
Weight ^(With int.bags)	5 ± 0,2 kg
Weight ^(Without int.bags)	3,2 ± 0,2 kg
Volume	45 l
Materials	Cordura, Al, Nylon, PE , TPU, Polyester
Loading capacity	30 kg

R-AID PRO - MULTIP. BACKPACK RED/BLACK W/5 POUCHES
CB04500A

Rev.3 (14/11/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it