

## R-AID PRO

### Multipurpose backpack



Multipurpose emergency backpack extremely versatile and robust. It is characterized by a structure integrated with the carrying points giving added strength and for a better support of the weight.

#### Specific features

Cordura fabric
Carrying point fixed to the supporting structure, with the internal part made of metal alloy and rigid polyethylene
TPU reinforcement on the bottom
2 external pockets with zip closure
1 removable top pocket with zip closure and equipped with elastic rings and pockets
1 removable front backpack attached by zip equipped with various straps and elastic rings, two shoulder straps and transparent front pocket
The padded main compartment provides maximum protection to the content and allows to store a set of bags
Strap bands to fix pockets or bags
Belts to fix an oxygen tank of 2 L
Ergonomic handles with loop for hanging kit and adjustable upper carrying point
Padded shoulder straps adjustable also in height in order to fit each operator
The padded lumbar support has a protection for the shoulder straps in order to protect them from dirt or damage
Anatomic padded back to provide maximum comfort
Reflective inserts increase visibility

#### Technical data

Dimensions	600 x 420 x 300 ± 30 mm
Weight	3,2 ± 0,2 kg
Volume	45 L
Materials	Cordura, Al, Nylon, PE, TPU, Polyester
Loading capacity	30 kg

R-AID PRO - MULTIP. BACKPACK RED/BLACK EMPTY  
**CB04550B**

Rev.4 (11/04/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)