

R-AID PRO

Multipurpose backpack



Multipurpose emergency backpack extremely versatile and robust. It is characterized by a structure integrated with the carrying points giving added strength and for a better support of the weight.

Specific features

- Cordura fabric
- Carrying point fixed to the supporting structure , with the internal part made of metal alloy and rigid polyethylene
- TPE reinforcement on the bottom
- 2 external pockets with zip closure
- 1 removable top pocket with zip closure and equipped with elastic rings and pockets
- 1 removable front backpack attached by zip equipped with various straps and elastic rings, two shoulder straps and transparent front pocket
- The padded main compartment provides maximum protection to the content and allows to store a set of bags
- Strap bands to fix pockets or bags
- Belts to fix an oxygen tank of 2lt
- Ergonomic handles with loop for hanging kit and adjustable upper carrying point
- Padded shoulder straps adjustable also in height in order to fit each operator
- The padded lumbar support has a protection for the shoulder straps in order to protect them from dirt or damage
- Anatomic padded back to provide maximum comfort
- Reflective inserts increase visibility

Technical data

| | |
|------------------|---|
| Dimensions | 600x420x300 ± 30 mm |
| Weight | 3,2 ± 0,2 kg |
| Volume | 45 l |
| Materials | Cordura, Al, Nylon, PE , TPU, Polyester |
| Loading capacity | 30 kg |

R-AID PRO - MULTIP. BACKPACK RED/BLACK EMPTY
CB04550A

Rev.1 (14/11/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it