

ALERT KRONE

Resuscitation case



Hard shell case for housing and transport devices generally used in resuscitation procedures.

Specific Characteristics

Extremely robust metal structure
Wide handle for easy transport
Metal support feet
Rounded corner protectors made of chromed steel
Rapid closure system with lock
180° opening to maximize access to content
Mobile internal dividers
The transparent separator allows easy view of the content and to keep it ordered
2 strap belt to fix an oxygen tank (max 1L)
3 rings for the attachment of the supplied straps allow to carry the case as a backpack

Technical data

External dimensions	525x165x375 ± 5 mm
Internal dimensions	495x150x355 ± 5 mm
Materials	AL, Acciaio, PE
Weight	5,7 ± 0,1 kg

ALERT KRONE - EMPTY-REINF.
ALUMIN. CASE
CB06001A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it