

LIENZ KRONE

Professional aluminum case
for advanced resuscitation equipment



Hard shell case for housing and transport devices generally used in resuscitation procedures.

Specific Characteristics

| |
|--|
| Extremely robust metal structure |
| Wide handle for easy transport |
| Metal support feet |
| Rounded corner protectors made of chromed steel |
| Rapid closure system with lock |
| 180° opening to maximize access to content |
| Mobile internal divisors |
| Thermo-formed compartments to store accessories, one of which is foldable |
| Dedicated compartments for oxygen tank and manual resuscitator |
| The transparent separator is fixed by magnet, allows easy view of the content and keeps it ordered |
| 2 strap belt to fix an oxygen tank (max 2 L without ring nut) |
| 3 rings for the attachment of the supplied straps allow to carry the case as a backpack |
| Equipped with phial holder |

Technical data

| | |
|---------------------|--------------------|
| External dimensions | 530x385x230 ± 5 mm |
| Materials | AL, Steel, PE |
| Weight | 6,5 ± 0,1 kg |

LIENZ KRONE - ALUMINIUM CASE
EMPTY
CB06002A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it