

ALERT KRONE

Resuscitation case



Hard shell case equipped with most commonly used devices during pediatric resuscitation procedures.

Specific Characteristics

Extremely robust metal structure
 Wide handle for easy transport
 Metal support feet
 Rounded corner protectors made of chromed steel
 Rapid closure system with lock
 180° opening to maximize access to content
 The transparent separator allows easy view of the content and to keep it ordered
 2 strap belt to fix an oxygen tank (max 1L)
 3 rings for the attachment of the supplied straps allow to carry the case as a backpack

Technical data

External dimensions	525x165x375 ± 5 mm
Internal dimensions	495x150x355 ± 5 mm
Materials	Al, Steel, PE
Empty case weight	5,7 ± 0,1 kg
Press.red. Output pressure	3,5/4 bar
Adjustable oxygen flow	0-1-1,5-2-3-4-6-9-12-15 l/min
Connection standard	UNI
Press. Red. output	Hose ø7mm

Included equipment

BG02007A	PRO B-LIFE- ADULT PVC MANUAL RESUSCITATOR
BS04011D	BURN S (BA) 10 X 10 CM - BURN STOP DRESSING
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS
DG01410C	DGX410 ADULT SPHYGMOMANOMETER
DG02110B	PHONENDOSCOPE
DG03100A	L-HANDLE CONVENTIONAL LIGHT
DG05001A	C-BLADE - SIZE 1-MACINTOSH BLADE
DG05002A	C-BLADE - SIZE 2-MACINTOSH BLADE
DG05003A	C-BLADE - SIZE 3-MACINTOSH BLADE
DG08062D	DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES
FG00607A	FORCEPS MAGILL CHILDREN
IF00150C	RUBBER TOURNIQUETTE
OX10005E	OXYGEN TANK 0,5 LT.
PR00308B	PRESS. RED. W/FLOW SELECT 15l/m
OX20120A	OXYMASK 50 MEDIUM CONCENTRATION
RM10321A	GUEDEL SIZE 60 MM – OROPGARYNGEAL CANN.
RM10322A	GUEDEL SIZE 70 MM – OROPGARYNGEAL CANN.
RM10323A	GUEDEL SIZE 80 MM – OROPGARYNGEAL CANN.
RM10892A	CUFFED MASK SZ.3 PVC COL. CODED WITHOUT VALVE
RM10893A	CUFFED MASK SZ 4
RM30400A	YOUNG - TONGUE FORCEPS
RV71101A	RES-Q-VAC ORANGE - MANUAL SUCTION PUMP
ST00630A	HT - ADULT ISOTHERMAL SHEET 220x160 cm
FG00605A	FORCEPS MAGILL ADULT

Med.Dev. CE marked according to the dir. 93/42/CE

ALERT KRONE COMPLETE ADULT-
 REINF. ALUMIN. CASE FOR RESUSC.
 EQUIP
CB06016A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it