

ALERT KRONE

Resuscitation case



Hard shell case equipped with most commonly used devices during pediatric resuscitation procedures.

Specific Characteristics

Extremely robust metal structure
Wide handle for easy transport
Metal support feet
Rounded corner protectors made of chromed steel
Rapid closure system with lock
180° opening to maximize access to content
The transparent separator allows easy view of the content and to keep it ordered
2 strap belt to fix an oxygen tank (max 1L)
3 rings for the attachment of the supplied straps allow to carry the case as a backpack

Technical data

External dimensions	525x165x375 ± 5 mm
Internal dimensions	495x150x355 ± 5 mm
Materials	AL, Acciaio, PE
Empty case weight	5,7 ± 0,1 kg
Press.red. Output pressure	3,5/4 bar
Adjustable oxygen flow	0-1-1,5-2-3-4-6-9-12-15 l/min
Connection standard	UNI
Press. Red. output	Hose ø7mm

Included equipment

BG02043A	PRO B-LIFE- INSUFF.PVC.PÉD. C/SURP.ET CUFFED MASK
BS04011D	BURN S (BA) 10 X 10 CM - BURN STOP DRESSING
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS
DG01414C	DGX414 W/ PAED. CUFF - ANEROID SPHYGM.
DG02110B	NURSE SCOPE - PHONENDOSCOPE
DG03101A	L-HANDLE CONVENTIONAL LIGHT PAEDIATRIC
DG05000A	C-BLADE - SIZE 0-MACINTOSH BLADE
DG05001A	C-BLADE - SIZE 1-MACINTOSH BLADE
DG05002A	C-BLADE - SIZE 2-MACINTOSH BLADE
DG08062D	DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES
FG00607A	FORCEPS MAGILL CHILDREN
IF00150C	RUBBER TOURNIQUETTE
OX10005E	OXYGEN TANK 0,5 LT.
PR00308B	PRESS. RED. W/FLOW SELECT 15l/m
OX20103A	OXYMASK PEDIATRIC MEDIUM CONCENTRATION
RM10325A	GUEDEL SIZE 40 MM – OROPGARYNGEAL CANN.
RM10320A	GUEDEL SIZE 50 MM – OROPGARYNGEAL CANN.
RM10321A	GUEDEL SIZE 60 MM – OROPGARYNGEAL CANN.
RM10892A	CUFFED MASK SZ.3 PVC COL. CODED WITHOUT VALVE
RM30300A	COLLIN - STAINLESS STEEL TONGUE FORCEPS
RV71101A	RES-Q-VAC ORANGE - MANUAL SUCTION PUMP
ST00621A	HT - PAED. ISOTHERMAL SHEET 140x90 cm

Med.Dev. CE marked according to the dir. 93/42/CE

ALERT KRONE COMPLETE PAEDIATRIC-
REINF. ALUMIN. CASE FOR RESUSC.
EQUIP.
CB06017A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it