

LIENZ KRONE

Professional aluminum case
for advanced resuscitation equipment



The image is only indicative



Hard shell case equipped with most commonly used devices during pediatric resuscitation procedures.

Specific Characteristics

- Extremely robust metal structure
- Wide handle for easy transport
- Metal support feet
- Rounded corner protectors made of chromed steel
- Rapid closure system with lock
- 180° opening to maximize access to content
- Thermo-formed compartments to store accessories, one of which is foldable
- Dedicated compartments for oxygen tank and manual resuscitator
- The transparent separator is fixed by magnet, allows easy view of the content and keeps it ordered
- 3 rings for the attachment of the supplied straps allow to carry the case as a backpack

Technical data

External dimensions	530x385x230 ± 5 mm
Materials	AL, Steel, PE
Empty case weight	6,5 ± 0,1 kg
Press.red. Output pressure	3,5/4 bar
Adjustable oxygen flow	0-1-1,5-2-3-4-6-9-12-15 l/min
Connection standard	UNI
Press. Red. output	Hose ø7mm

Included equipment

BG02043A	PRO B-LIFE- INSUFF.PVC.PÉD. C/SURP.ET CUFFED MASK
BS04011D	BURN S (BA) 10 X 10 CM - BURN STOP DRESSING
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS
DG01414C	DGX414 W/ PAED. CUFF - ANEROID SPHYGM.
DG02110B	NURSE SCOPE - PHONENDOSCOPE
DG03101A	L-HANDLE CONVENTIONAL LIGHT PAEDIATRIC
DG05000A	C-BLADE - SIZE 0-MACINTOSH BLADE
DG05001A	C-BLADE - SIZE 1-MACINTOSH BLADE
DG05002A	C-BLADE - SIZE 2-MACINTOSH BLADE
DG08062D	DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES
FG00607A	FORCEPS MAGILL CHILDREN
IF00150C	RUBBER TOURNIQUETTE
OX10039E	OXYGEN TANK 2 LT.
PR00308B	PRESS. RED. W/FLOW SELECT 15l/m
OX20103A	OXYMASK PEDIATRIC MEDIUM CONCENTRATION
RM10325A	GUEDEL SIZE 40 MM – OROPGARYNGEAL CANN.
RM10320A	GUEDEL SIZE 50 MM – OROPGARYNGEAL CANN.
RM10321A	GUEDEL SIZE 60 MM – OROPGARYNGEAL CANN.
RM10892A	CUFFED MASK SZ.3 PVC COL. CODED WITHOUT VALVE
RM30300A	COLLIN - STAINLESS STEEL TONGUE FORCEPS
RV71101A	RES-Q-VAC ORANGE - MANUAL SUCTION PUMP
ST00621A	HT - PAED. ISOTHERMAL SHEET 140x90 cm

Med.Dev. CE marked according to the dir. 93/42/CE

LIENZ KRONE COMPLETE PAEDIATRIC-
PROFESS. ALUM. CASE FOR ADV.
RESUSC. EQUIP.
CB06019A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it