

# KDX 1

## Resuscitation case



Hard shell case for housing and transport of devices generally used in resuscitation procedures.

### Specific Characteristics

- Quick release closure system
- Manual pressurizin valve
- Extremely robust structure
- Seal to provide protection against ingress of fluids or dirt
- Internal foam for maximum protection to the content
- Thanks to the pre-made cutting of the foam, the user can customize the housing as needed.
- Predisposed for fixation of metal supports

### Technical data

External dimensions	405x515x200 ± 5 mm
Internal dimensions	355x475x180 ± 5 mm
Materials	PP, PE
Weight without foam	3,6 ± 0,1 kg
Wight with foam	4 ± 0,1 kg

KDX 1 - KDX 1 EMPTY -  
RESUSCITATION CASE  
**CB06045C**

Rev.2 (24/08/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)