

# PRO KDX

## Resuscitation case



Hard shell case equipped with devices generally used in resuscitation procedures.

### Specific Features

Quick release closure system  
 Manual pressurizing valve  
 Extremely robust structure  
 Sealed to provide protection against ingress of fluids or dirt

### Technical data

External dimensions	405x515x200 ± 5 mm
Internal dimensions	355x475x180 ± 5 mm
Materials	PP, PE

### Included equipment

DG03100 – ADULT LARYNGOSCOPE HANDLE  
 DG05001 – C-BLADE MACINTOSH BLADE SZ 1  
 DG05002 – C-BLADE MACINTOSH BLADE SZ 2  
 DG05003 – C-BLADE MACINTOSH BLADE SZ  
 BG02041 – PRO B-LIFE ADULT WITH VALVE AND MASK SZ 4  
 RM10892 – CUFFED MASK SZ 3  
 FG00607 – MAGILL PAEDIATRIC - TUBE-GUIDING FORCEPS  
 OX20120 – OXYMASK 50 ADULT MEDIUM CONC. OXYG.MASK  
 DG02110 – DGX 110 – NURSE SCOPE - PHONENDOSCOPE  
 IF00150 – IFX 50 RUBBER TOURNIQUETTE  
 ST00630 - HT ADULT - ISOTHERMAL SHEET  
 FG00605 - MAGILL ADULT - TUBE-GUIDING FORCEPS  
 BS04011 - BURN S (BA) 10 X 10 CM - BURN STOP DRESSING  
 RM30400 - YOUNG - STAINLESS STEEL TONGUE FORCEPS  
 PR00308 - PRESSURE REDUCER W/FLOW SELECT 15l/m  
 OX10020 – UNI OXYGEN TANK 2 LT.  
 RV71150 - B-VAC - MANUAL ASPIRATOR FOR SECRETIONS  
 CB02030 – T2 BAG  
 DG00603 - UNIVERSAL BLACK - TUFF CUT SCISSORS  
 RM10321 - GUEDEL - DISPOSABLE OROPH. AIRWAY 60 MM  
 RM10322 - GUEDEL - DISPOSABLE OROPH. AIRWAY 70 MM  
 RM10323 - GUEDEL - DISPOSABLE OROPH. AIRWAY 80 MM  
 DG08062 - DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES  
 DG01410 - DGX 410 W/ AD. CUFF - ANEROID SPHYGMOMANOMETER

Med. Dev. CE marked according to dir. 93/42/CEE

PRO KDX - WATERTIGHT  
 RESUSCITATION CASE  
**CB06073C**



Rev.1 (03/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)