

PRO KDX

Resuscitation case



Hard shell case equipped with devices generally used in resuscitation procedures.

Specific Features

Quick release closure system
 Manual pressurizing valve
 Extremely robust structure
 Sealed to provide protection against ingress of fluids or dirt

Technical data

External dimensions	405x515x200 ± 5 mm
Internal dimensions	355x475x180 ± 5 mm
Materials	PP, PE

Included equipment

DG03100 – ADULT LARYNGOSCOPE HANDLE
 DG05001 – C-BLADE MACINTOSH BLADE SZ 1
 DG05002 – C-BLADE MACINTOSH BLADE SZ 2
 DG05003 – C-BLADE MACINTOSH BLADE SZ
 BG02041 – PRO B-LIFE ADULT WITH VALVE AND MASK SZ 4
 RM10892 – CUFFED MASK SZ 3
 FG00607 – MAGILL PAEDIATRIC - TUBE-GUIDING FORCEPS
 OX20120 – OXYMASK 50 ADULT MEDIUM CONC. OXYG.MASK
 DG02110 – DGX 110 – NURSE SCOPE - PHONENDOSCOPE
 IF00150 – IFX 50 RUBBER TOURNIQUETTE
 ST00630 - HT ADULT - ISOTHERMAL SHEET
 FG00605 - MAGILL ADULT - TUBE-GUIDING FORCEPS
 BS04011 - BURN S (BA) 10 X 10 CM - BURN STOP DRESSING
 RM30400 - YOUNG - STAINLESS STEEL TONGUE FORCEPS
 PR00308 - PRESSURE REDUCER W/FLOW SELECT 15l/m
 OX10020 – UNI OXYGEN TANK 2 LT.
 RV71150 - B-VAC - MANUAL ASPIRATOR FOR SECRETIONS
 CB02030 – T2 BAG
 DG00603 - UNIVERSAL BLACK - TUFF CUT SCISSORS
 RM10321 - GUEDEL - DISPOSABLE OROPH. AIRWAY 60 MM
 RM10322 - GUEDEL - DISPOSABLE OROPH. AIRWAY 70 MM
 RM10323 - GUEDEL - DISPOSABLE OROPH. AIRWAY 80 MM
 DG08062 - DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES
 DG01410 - DGX 410 W/ AD. CUFF - ANEROID SPHYGMOMANOMETER

Med. Dev. CE marked according to dir. 93/42/CEE

PRO KDX - WATERTIGHT
 RESUSCITATION CASE
CB06073C



Rev.1 (03/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it