


## RESPONSE

### Multi purpose emergency bag



Professional bag ideal to store rescue equipment and accessories, characterized by easy access to the content thanks to the independent compartments.

#### Specific features

Rigid container in the lower part dedicated to fix an oxygen tank by means of 3 straps and able to store accessories

Upper compartment with access from the top through a flap with double zip closure

2 elastic bands to store accessories on the internal side of the closure flap

1 elastic band and large pocket with strap closure on the internal side of the upper compartment

Transparent window on the outside with access from the inside

Shoulder strap and handles allow easy transport

Coated bottom on the outside

#### Included equipment

2 Matrix bags – Yellow phial holder, orange dedicated to infusion equipment

#### Technical data

Dimensions	480x260x320 ± 10 mm
Volume	23 ± 0,5 l
Shoulder strap length	Up to 100 cm
Bag weight	1,99 ± 0,1 kg
Weight 1 Matrix bag	380 g
Materials	Nylon, PVC, PE, EVA

PRO RESPONSE – PROF EMERGENCY BAG  
W/2 MATRIX BAGS  
CB07021A

Rev.0 (05/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)