

## LIFE BAG 2

### Rescue bag



Professional bag ideal to store first aid devices and accessories, characterized by quick and easy access to the content thanks to wide openings and independent compartments.

#### Specific features

Equipped with handles and detachable shoulder strap  
 Main compartment with 3-compartment separator and top flap with elastic bands to attach accessories  
 The internal separator is detachable and has an integrated carrying handle  
 2 side compartments and a front one with zip closure  
 The padding guarantees good content protection  
 PVC reinforced bottom

#### Technical data

Dimensions	420x240x200 ± 20mm
Overall volume	11,5 ± 0,3 l
Weight	890 ± 10 g
Materials	Nylon, PVC

LIFE BAG 2 – RESCUE BAG  
 CB07040A

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)