

LIFE BAG 2

Rescue bag



Professional bag ideal to store first aid devices and accessories, characterized by quick and easy access to the content thanks to wide openings and independent compartments.

Specific features

- Equipped with handles and detachable shoulder strap
- Main compartment with 3-compartment separator and top flap with elastic bands to attach accessories
- The internal separator is detachable and has an integrated carrying handle
- 2 side compartments and a front one with zip closure
- The padding guarantees good content protection
- PVC reinforced bottom

Technical data

Dimensions	420x240x200 ± 20mm
Overall volume	11,5 ± 0,3 l
Weight	890 ± 10 g
Materials	Nylon, PVC

LIFE BAG 2 – RESCUE BAG
CB07040A

Rev.1 (15/10/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it