

Emergency Emergency bag



Practical bag with everything is needed for an efficient approach to primary resuscitation.

Specific Characteristics

- High portability thanks to the shoulder strap
- High visibility due to reflective bands
- Compact and practical to use
- Logical and ordered placing of the content thanks to the rigid separation of the main compartment
- 2 lateral pockets with zip closure
- Reinforced bottom

Equipment

FG00605A	Magill adult forceps
DG02110B	Black phonendoscope
ST00630A	HT Isothermic sheet adult
IF00150C	IFX50 Tubular tourniquet
OX20130A	OXX30 O ₂ anticrush tube
RM10892A	Cuffed mask mis.3
RM10893A	Cuffed mask mis.4
RM30700A	Moulder jaw spreader
OX10005E	O ₂ bottle 0,5l
PR00304B	Pressure reducer
RM30400A	Young steel tongue holder
DG00603A	Tuff cut Universal scissors
DG01411C	Sphygmomanometer
RM10321A	Disposable non-sterile Guedel sz.1
RM10322A	Disposable non-sterile Guedel sz.2
RM10323A	Disposable non-sterile Guedel sz.3
RM10324A	Disposable non-sterile Guedel sz.4
RV71120B	RES-Q-VAC Manual suction device
BG02007A	Pro-B-Life adult PVC disposable resuscitator
DG08062D	DGX062 Pen Light blue no/batt

Technical data

Dimensions	490x230x220 mm
Bag material	Nylon
Bag weight	730 g
Kit weight	4,8 kg
Tank connection standard	UNI
Reducer outgoing pressure	3,5 bar
Reducer outgoing flow	6 l/min
Anticrush tube length	2,1 m
Suction device max vacuum	300 mmHg
Suction device Max flow	22 l/min
Manual resuscitator volume	2500 ml
Man. resusc. Overp. Valve calibration	60 mbar

EMERGENCY – COMPLETE
EMERGENCY BAG
CB08000A



Rev.1 (05/11/2014)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it