

Emergency Emergency bag



Practical bag ideal to store and carry the devices used in primary resuscitation procedures. Available in many other languages.

Specific Characteristics

- High portability thanks to the shoulder strap
- High visibility due to reflective bands
- Compact and practical to use
- Rigid separation of the main compartment
- 2 lateral pockets with zip closure
- Reinforced and coated bottom makes the cleaning easier
- Possibility to store an oxygen tank (max 1l without protection)

Technical data

Dimensions	490x230x220 mm
Bag material	Nylon, PE
Bag weight	700 ± 30 g

EMERGENCY – EMPTY BAG WITH
SPANISH WRITTEN
CB08006A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it