

## BLUE BAG 4

### Backpack/bag



Blue bag 4 is a large backpack, also usable as a bag, designed for transport of an oxygen tank.

#### Specific features

The front compartment can be completely opened by a zip closure and has 2 detachable pockets, 2 elastic bands and 3 pockets with strap closure to store accessories

The main compartment is padded and has a zip closure. The paddings on the back side have a shape designed to store an oxygen tank up to 3 liters of capacity. The oxygen tank is fastened by means of two bands with adjustable strap

Side pocket with zip closure

Reflective inserts to increase visibility

Padded backrest to increase comfort of the rescuer

In the back there is a PVC compartment where is possible to place the shoulder straps when not in use

Carrying handle to improve portability

#### Technical data

Maximum dimensions	560x360x260 ± 10mm
Overall volume	43 ± 0,5 l
Internal bags	2
Bags color	Yellow, Red
Materials	Nylon, PVC
Weight	1200 ± 20g
Loading capacity	15 kg

BLUE BAG 4 - BACKPACK/BAG

CB09001C

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)