

OXYPACK 2

Oxygen tank carrier



Oxygen tank carrier design for transport of 2 liters oxygen tanks.

Specific features

- Easy to clean PVC sheet
- Adjustable shoulder strap to allow easy transport
- Hole between the two compartments to allow access to a flowmeter that could be connected to the oxygen tank
- Front pocket with zip closure to store accessories used for oxygen administration
- Padded main compartment to protect the content and to allow a more comfortable transportation

Technical data

| | |
|-----------------------|-----------------------|
| Dimensions | 150x290x500 ± 10 mm |
| Pocket dimensions | 370x100x130 ± 10 mm |
| Shoulder strap length | From 60 to 120 ± 2 cm |
| Loading capacity | 6 Kg |
| Weight | 490 ± 10 g |

OXYPACK 2 - PVC OXYGEN
TANK FOR TANK 2 L
CB09002A

Rev.3 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it