

# OXYPACK 1

## Oxygen tank carrier



Oxygen tank carrier design for transport of 5 liters oxygen tanks.

### Specific features

- Easy to clean PVC sheet
- Adjustable shoulder strap to allow easy transport
- Hole between the two compartments to allow access to a flowmeter that could be connected to the oxygen tank
- Front pocket with zip closure to store accessories used for oxygen administration
- Padded main compartment to protect the content and to allow a more comfortable transportation

### Technical data

Tank compartment diameter	150 ± 10 mm
Height	640 ± 10 mm
Dimensions lying bag <sup>Empty</sup>	650x275x80mm approx
Pocket dimensions	370x130x80 ± 10mm
Shoulder strap length	From 60 to 120 ± 5 cm
Weight	610 ± 10 g

OXYPACK BAG FOR 5 LT.  
 OXYGEN TANK EMPTY  
**CB09005A**

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)