

SAM

Bag for oxygen tank



Bag designed to carry an oxygen tank and other devices commonly used in ventilation procedures.

Specific features

- Handles and shoulder strap for easy transport
- One main compartment with zip closure and small window to allow access to the oxygen tank valve or to the attached flowmeter
- 2 lateral wide pockets with zip closure
- Can contain oxygen tanks up to 5 liters

Technical data

Dimensions	650x250x240 ± 20 mm
Material	PVC, Nylon
Weight empty	2,3 ± 0,1 kg

SAM - PVC OXYGEN TANK BAG EMPTY BLUE
CB09010A

Rev.1 (20/02/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it