

OXYPACK S

Oxygen tank carrier for stretcher



Oxygen tank carrier designed for transport of 2 liters oxygen tanks to be attached to the stretcher.

Specific features

Easy to clean PVC sheet

Adjustable shoulder strap to allow easy transport

Hole between the two compartments to allow access to a flowmeter that could be connected to the oxygen tank

Front pocket with zip closure to store accessories used for oxygen administration

Padded main compartment to protect the content and to allow a more comfortable transportation

Quick release buckles for attachment to the stretcher

Technical data

Dimensions	150x290x500 ± 10 mm
Pocket dimensions	370x100x130 ± 10 mm
Shoulder strap length	From 60 to 120 ± 2 cm
Loading capacity	6 Kg
Weight	550 ± 10 g

**OXYPACK S – OTXYGEN TANK
CARRIER FOR STRETCHERS
CB09012A**

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it