

OXYPACK S

Oxygen tank carrier for stretcher



Oxygen tank carrier designed for transport of 2 liters oxygen tanks to be attached to the stretcher.

Specific features

- Easy to clean PVC sheet
- Adjustable shoulder strap to allow easy transport
- Hole between the two compartments to allow access to a flowmeter that could be connected to the oxygen tank
- Front pocket with zip closure to store accessories used for oxygen administration
- Padded main compartment to protect the content and to allow a more comfortable transportation
- Quick release buckles for attachment to the stretcher

Technical data

Dimensions	150x290x500 ± 10 mm
Pocket dimensions	370x100x130 ± 10 mm
Shoulder strap length	From 60 to 120 ± 2 cm
Loading capacity	6 Kg
Weight	550 ± 10 g

**OXYPACK S – OTXYGEN TANK
CARRIER FOR STRETCHERS
CB09012A**

Rev.0 (10/02/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it