

## TANKER

Oxygen tank carrier for stretcher side bar



Oxygen tank carrier with shoulder strap designed for mounting on the side bar of Spencer stretchers.

### Specific features

- Robust polyethylene body easy to clean
- Adjustable shoulder strap for an easy transport
- The device can be anchored to the side of the stretcher using the dedicated support sold separately
- The holes on the body allow the view and adjustment of oxygen tank controls (manometer, valve, flowmeter)
- Designed to contain a 2liters oxygen tank with reducing valve

### Optional accessories

CB09028C- TANKER SUPPORT FOR SAFEBAR

### Technical data

Dimensions	430x180x160 ± 5 mm
Depth	390mm <small>Space for valve included</small>
Internal diameter	115 mm
Materials	PE, Nylon
Weight	950 ± 10g

Class I med.dev. according to Dir. 93/42/CEE

TANKER YELLOW OXYGEN  
TANK HOLDER FOR  
STRETCHER'S SIDEBARS  
FIXATION

CB09025C

CND Classification



V9099

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)