

TANKER SUPPORT

For Safe Bar



Support for oxygen tank holder Tanker to be installed on sidebars Safe Bar.

Specific features

- Robust polyethylene structure easy to clean
- Easy to install
- The tank holder is easily inserted by means of its feet into the dedicated holes and acting on the locking system with push-button disengagement

Technical data

Dimensions	350x165x100 ± 5 mm
Materials	PE, Nylon, Steel
Weight	580 ± 10g

Class I Med. Dev. according to Dir. 93/42/CEE

TANKER SUPPORT FOR
SIDEBAR

CB09028C

CND Classification



V9099

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it