

OXYPACK S

Oxygen tank carrier for stretcher



Oxygen tank carrier designed for transport of 5 liters oxygen tanks to be attached to the stretcher.

Specific features

- Easy to clean PVC sheet
- Adjustable shoulder strap to allow easy transport
- Hole between the two compartments to allow access to a flowmeter that could be connected to the oxygen tank
- Front pocket with zip closure to store accessories used for oxygen administration
- Padded main compartment to protect the content and to allow a more comfortable transportation
- Quick release buckles for attachment to the stretcher

Technical data

Tank compartment diameter	150 ± 10 mm
Dimensions lying bag ^{Empty}	650x275x80mm approx
Pocket dimensions	370x130x80 ± 10mm
Shoulder strap length	From 60 to 120 ± 2 cm
Loading capacity	6 Kg
Weight	760 ± 10 g

**OXYPACK S – OXYGEN TANK
CARRIER FOR STRETCHERS
CB09036A**

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it