BLUE BAG 2 Professional bag





Professional bag particularly suitable to store the main devices commonly used for measurement of vital signs.

Specific features

Frontal pocket with zip closure

- Two side pockets with zip closure
- Top pocket with zip closure
- Main compartment with zip closure and equipped with:
- 5 internal bags of different colors with zip closure and transparent window
- Pocket with zip closure on the internal side of the frontal flap
- Internal separator characterized by a side made of tissue for the
- fixing of the internal bags; on the opposite side there are $\ensuremath{\mathsf{3}}$
- transparent compartments with zip closure to store additional accessories
- The main compartment is padded to provide protection to the content
- Pocket on the back with strap closure
- Maximum portability thanks to the handles and detachable shoulder strap

Reflective inserts to increase visibility

Technical data	
Dimensions	390x210x340 ± 10mm
Volume	20 ± 0,5 l
Shoulder strap lenght	From 65 to 100cm
N° of bags	5
Bags color	Green, orange, red, yellow, blue
Materials	Nylon, PVC
Weight	$1020 \pm 20g$

BLUE BAG 2 - PROFESSIONAL BAG CB09601C

Rev.1 (15/01/2018) UNCHECKED COPY – further revisions will be available on <u>http://support.spencer.it</u>

Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222