

BLUE BAG 2

Complete professional bag



Professional bag equipped with devices commonly used in first aid procedures.

Specific features

Frontal pocket with zip closure
Two lateral pockets with zip closure
Top pocket with zip closure
Main compartment with zip closure and equipped with
- 5 internal color coded pockets with zip closure and transparent top side
- Pocket with zip closure placed on the internal side of the front flap
- Internal separator characterized by a fabric side to fix the pockets and another one with 3 transparent compartments with zip closure. Between the two sides there is a compartment to store additional accessories
Padded main compartment for maximum protection to the content
Pocket on the back side with zip closure
Maximum portability thanks to the carrying handles and the detachable shoulder strap
Reflective bands increase visibility

Technical data

Dimensions	390x210x340 ± 10mm
Volume	20 ± 0,5 l
Shoulder strap length	From 65 to 100cm
Internal bags	5
Internal bags colours	Green, orange, red, yellow, blue
Materials	Nylon, PVC
Empty bag weight	980 ± 20g

Equipment

BS04011D	2x BURN S (BA) 10 X 10 CM - BURN STOP DRESSING
BS04021D	2x BURN S (BA) 20 X 20 CM - BURN STOP DRESSING
DG00602A	UNIVERSAL YELLOW - TUFF CUT SCISSORS
DG01418A	PROFESSIONAL ANEROID SPHYGMOMANOMETER (AD.)
DG08062D	DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES
FG03384B	KOCHER HEMOSTATIC FORCEPS
GM00747K	4x TRIANGULAR BAND.FOR CLAVICL E FRACT. TNT CM.9
GM02101K	4x GAUZE COMPR.STER. 10X10CM 25 SHEETS/pack
GM02106K	2x STERILE OCULAR GAUZE COMPRESS
GM02250K	2x CLOTH PLASTER 5CM X 5M
GM02405K	2x GAUZE BENDAGE 5M x5CM
GM02407K	2x GAUZE BANDAGE 5M X 7CM
GM02410K	2x GAUZE BANDAGE 5M X 10CM
GM02607K	4x ELASTIC BANDAGE FOR MEDICATION 7 cm x 4,5 m
GM02610K	4x ELASTIC BANDAGE FOR BANDAGING 10 cm x 4,5 m
GM04250K	GMX003 – DISINFECTANT 250ML
IF00150C	RUBBER TOURNIQUETTE
IF00426A	2x SPENCER ICE PACK PE
RM10322A	GUEDEL - DISPOSABLE OROPH. AIRWAY 70 MM
RM10323A	GUEDEL - DISPOSABLE OROPH. AIRWAY 80 MM
RM10324A	GUEDEL - DISPOSABLE OROPH. AIRWAY 90 MM
RM10327A	GUEDEL - DISPOSABLE OROPH. AIRWAY 100 MM
ST00630A	ISOTHERMIC SHEET ADULT SIZES 160X220 CM.
TA09019A	FLEX MASK SET WITH ORANGE PVC CASE

Med. Dev. CE marked according to dir.93/42/CEE

MOD.E020 BLUE BAG 2 COMPLETE
CB09609

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it