


Command

Special backpack for fire fighters



Backpack designed for use by fire fighters, ideal for storage and transport of a lifting harness or an Oxygen tank.

Specific features

- Lifting handle fixed on aluminum reinforcement
- Its regular shape allows to optimize the bulk maximizing the volume
- 2 side pockets with zip closure. Inside these pockets, strap surfaces allow fixation of R-Series bags.
- Main compartment with 4 side mesh pockets and belts to fix the oxygen tank
- On the front flap, 2 pockets with zip closure allow to store accessories
- Adjustable and padded shoulder straps
- Adjustable lumbar belt for a better distribution of the load
- Backrest with aluminum reinforcement
- Reflective inserts to increase visibility
- PVC reinforced bottom to avoid permeability of the backpack when placed on the ground

Standard equipment

R-Series Bags	
CB90003	Sacca R1 – Suitable for Intubation
CB90005	Sacca R2 – Suitable for ventilation
CB90002	Sacca R3 – Suitable for phial holder
CB90004	Sacca R4 – Suitable for infusion

Technical data

Dimensions	300x400x530 mm
Weight without bags	2,45 ± 0,1 kg
Weight with bags	3,50 ± 0,2 kg
Volume	44,5 l
Materials	Nylon, PVC, Al
Loading capacity	20 kg
Oxygen tank max. dimensions	ø220x470 ± 30 mm
Can store Spencer oxygen tanks up to 2liters	

COMMAND - BACKPACK RED SPEC.FIRE FIGHT.W/4
POUCH.
CB90001C

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it