

## R SERIES

### Phial holder bag R3



Bag designed for use with Spencer backpacks or bags equipped with fixing straps. It can also be used separately and is ideal for use as phial holder.

#### Specific Characteristics

Strap band on rear side for a fixed storage inside Spencer backpack or bags

Transparent window on the front side to place useful information about the content

Carrying handle

Two internal transparent pockets with zip closure allow excellent space allocation and allow to view the content

Padded divisor for better protection during transport

Equipped with adjustable straps to ensure a firm holding of phials according to their number and size.

#### Technical data

Dimensions	220x230x95 mm
Weight	270 g
Materials	Nylon, PVC
Storable phials	22 small in the adjustable straps 14 medium in the opposite side

R 3 – ORANGE PHIAL HOLDER BAG

CB90002A

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)