

CROSS UP 8409

Stretcher with telescopic handles



Mattress shown in the image is sold separately

Cross Up 8409 is a stretcher for recovery and transport of patients to be used alone or with the dedicated trolley. The trolley allows a secure fixing inside a medical vehicle.

Specific Characteristics

Polyethylene board easy to clean
Adjustable backrest
Telescopic handles
Can be used with the Cross Up 8409 trolley
Light and robust
Foldable sidebars <i>Safe bar</i>
Wheels \varnothing 100 mm to increase sliding action
Trendelenburg antishock position

Standard equipment

ST70002A	Black belt 2 pcs, metal buckle, reflex
----------	--

Accessories

CR00007E	Trolley Cross Up 8409
CR00008E	Trolley Cross Up 8409 Twist
CR60003E	Trolley Cross Up 8409 Cert. 10G
CR60004E	Trolley Cross Up 8409 Twist Cert. 10G
ST70000A	Anatomic PVC mattress QMX/777
IF01043B	IV Holder TRACK 4
ST00497B	DNA strap belt w/integrated re-winding system
ST00498B	DNA strap thoracic belt w/re-winding system
CB09028C	Tanker fixation support for safe bar
ST00491A	STX 90 telescopic head rest for stretchers

Technical data

Width	560 \pm 10 mm
Length with handles closed	1970 \pm 10 mm
Length with handles opened	2500 \pm 10 mm
Backrest positions	5
Wheels	\varnothing 100 mm
Sidebars length	680 mm
Sidebars height	200 mm
Materials	Al, Steel, PE
Weight	21 kg
Loading capacity	250 kg

Class I Med.Device according to Dir. 93/42/CEE

CROSS UP STRETCHER 8409 T

CR00010E

CND Classification



V0899

Rev.0 (12/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it