

CROSS UP 8409

Stretcher with telescopic handles



The mattress shown in the picture is sold separately. The image could differ from the real product

Cross Up 8409 is a stretcher for recovery and transport of patients to be used alone or with the dedicated trolley which allows easy fixation inside the ambulance using Spencer fasteners.

Specific features

- Poliethylene board easy to clean
- Adjustable backrest
- Telescopic handles
- Can be used with the trolley Cross Up 8409
- Robust and lightweight
- Foldable sidebars *Safe bar*
- Ø100 mm wheels to improve handling
- Anti-shock positions Trendelenburg and Fowler

Standard equipment

ST70002A	Belt 2 pieces metal reflex black
----------	----------------------------------

Optional accessories

CR00007E	Cross Up Trolley 8409
CR00008E	Cross Up Trolley 8409 Twist
CR60003E	Cross Trolley Up 8409 Cert. 10G
CR60004E	Cross Up Trolley 8409 Twist Cert. 10G
ST70000A	Anatomic mattress QMX/777
IF01043B	IV Holder TRACK 4
ST00497B	DNA Strap w/integr. rewinding system
ST00498B	DNA Strap torax. w/int. Rewinding system
CB09028C	Tanker fixation support for safe bar
ST00491A	STX 90 telescopic head rest for stretchers

Technical data

Width	560 ± 10 mm
Length with handles closed	1970 ± 10 mm
Length with handles opened	2500 ± 10 mm
Backrest positions	5
Wheels	Ø100 mm
Sidebars length	680 mm
Sidebars height	200 mm
Materials	Al, Stainless steel, PE
Weight	21 kg
Loading capacity	250 kg

Class I Med.Device according to Dir. 93/42/CEE

STRETCHER CROSS UP 8409 TF

CR00012E

CND Classification



V0899

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it