

## CROSS UP 8409

### Stretcher with housing for spine board



Mattress and spinal board shown in the image are sold separately

Cross Up 8409 is a stretcher for recovery and transport of patients to be used alone or with the dedicated trolley. The trolley allows a secure fixing inside a medical vehicle. In this version, the stretcher is equipped with a housing for a B-Bak spinal board. This technical solution, allow to reduce considerably the time needed to carry out recovery operations.

#### Specific Characteristics

Dedicated housing for B-Bak spinal board  
 Polyethylene board easy to clean  
 Adjustable backrest  
 Telescopic handles  
 Can be used with the Cross Up 8409 trolley  
 Light and robust  
 Foldable sidebars *Safe bar*  
 Wheels  $\varnothing$ 100 mm to increase sliding action

#### Standard equipment

ST70002A Black belt 2 pcs, metal buckle, reflex

#### Accessories

CR00007E	Trolley Cross Up 8409
CR00008E	Trolley Cross Up 8409 Twist
CR60003E	Trolley Cross Up 8409 Cert. 10G
CR60004E	Trolley Cross Up 8409 Twist Cert. 10G
ST02061B	B-BAK YELLOW spinal board with pin
ST70000A	Anatomic PVC mattress QMX/777
IF01043B	IV Holder TRACK 4
ST00497B	DNA strap belt w/integrated re-winding system
ST00498B	DNA strap thoracic belt w/re-winding system
CB09028C	Tanker fixation support for safe bar
ST00491A	STX 90 telescopic head rest for stretchers

#### Technical data

Width	560 $\pm$ 10 mm
Length with handles closed	1970 $\pm$ 10 mm
Length with handles opened	2500 $\pm$ 10 mm
Backrest positions	5
Wheels	$\varnothing$ 100 mm
Sidebars length	680 mm
Sidebars height	200 mm
Materials	Al, Steel, PE
Weight	21 kg
Loading capacity	250 kg

Class I Med.Device according to Dir. 93/42/CEE

STRETCHER 8409 S FOR B-BAK PIN

CR00013E

CND Classification



V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)