

CROSS UP 8409

Stretcher with housing for spine board



Mattress and spinal board shown in the image are sold separately

Cross Up 8409 is a stretcher for recovery and transport of patients to be used alone or with the dedicated trolley. The trolley allows a secure fixing inside a medical vehicle. In this version, the stretcher is equipped with a housing for a B-Bak spinal board. This technical solution, allow to reduce considerably the time needed to carry out recovery operations.

Specific Characteristics

- Dedicated housing for B-Bak spinal board
- Polyethylene board easy to clean
- Adjustable backrest
- Telescopic handles
- Can be used with the Cross Up 8409 trolley
- Light and robust
- Foldable sidebars *Safe bar*
- Wheels $\varnothing 100$ mm to increase sliding action
- Trendelenburg and Fowler antishock positions

Standard equipment

ST70002A	Black belt 2 pcs, metal buckle, reflex
----------	--

Accessories

CR00007E	Trolley Cross Up 8409
CR00008E	Trolley Cross Up 8409 Twist
CR60003E	Trolley Cross Up 8409 Cert. 10G
CR60004E	Trolley Cross Up 8409 Twist Cert. 10G
ST02061B	B-BAK YELLOW spinal board with pin
ST70000A	Anatomic PVC mattress QMX/777
IF01043B	IV Holder TRACK 4
ST00497B	DNA strap belt w/integrated re-winding system
ST00498B	DNA strap thoracic belt w/re-winding system
CB09028C	Tanker fixation support for safe bar
ST00491A	STX 90 telescopic head rest for stretchers

Technical data

Width	560 \pm 10 mm
Length with handles closed	1970 \pm 10 mm
Length with handles opened	2500 \pm 10 mm
Backrest positions	5
Wheels	$\varnothing 100$ mm
Sidebars length	680 mm
Sidebars height	200 mm
Materials	Al, Steel, PE
Weight	21 kg
Loading capacity	250 kg

Class I Med.Device according to Dir. 93/42/CEE

STRETCHER 8409 TF FOR B-BAK PIN

CR00015E

CND Classification



V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it