

## CROSS UP 8409

### 10G Self-loading trolley with removable stretcher



Spencer Cross Up 8409 is a self-loading trolley with removable stretcher and adjustable heights for patient transport. It is usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically.

#### Specific features

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
Antishock position Trendelenburg
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
Adjustable heights system allows an easier placing of the patient on the stretcher
The stretcher is easily unlocked from the trolley and is equipped with telescopic handles to allow easier lifting
10G certified if used with fasteners S-MAX or R-MAX

#### Standard Equipment

ST70000A	PVC anatomic mattress QMX/777
ST70002A	STX 702 two pieces black reflex metal belts

#### Technical data

Length	1975 ± 10mm
Width	580 ± 10 mm
Trolley weight	31 ± 1 kg
Stretcher weight	21 ± 1 kg
Loading capacity	250 kg
Wheels	ø200 mm <sup>±5%</sup>
Height loading carriage wheel	65 cm
Materials	Stainless steel, Al, PE
Side bars length	680 mm
Side bars Height	200 mm

#### Accessories

ST42702C	Fastener R-Max <small>(Mandatory for 10G versions)</small>
ST42707B	Fastener S-Max <small>(Mandatory for 10G versions)</small>
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integr. re-winding system
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

#### Class I Med.Device according to Dir. 93/42/CEE

CROSS UP 8409 W/STRETCHER T  
W/QMX777 CERT. 10G  
**CR00022E**



<b>CND Classification</b>	V0899
<b>Registration number</b>	1110306

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)