

# CROSS CHAIR

## Stretcher/Chair



Cross Chair is a device that adopts innovative solutions for the transport of patients to be used in rescue operations on road and ambulances. The device can be used in chair or stretcher configuration, adapting to wider rescue scenarios. To reach more versatility and safety, there is a dedicated trolley to allow the use on ambulances and make transport easier.

### Specific Characteristics

- Configuration change takes less than 2 seconds
- 4-segment PVC Mattress with removable cover
- Polyethylene board easy to clean
- Backrest adjustable in 5 positions
- Telescopic handles
- Auto positioning rear handles
- Suitable for use with Cross Chair Trolley
- Light and robust
- The side bars *Safe bar* contain the patient also in Chair configuration
- Wheels  $\varnothing 150$  mm to increase sliding action
- Trendelenburg antishock position
- Swiveling front  $\varnothing 75$  mm wheels with brakes that lock rolling and axis rotation make transport in chair position easier

### Standard equipment

ST70002A	Black belt 2 pcs, metal buckle, reflex
ST43014	Cross Chair mattress

### Accessories

IF01043B	IV holder TRACK 4
CR00081	Cross Chair trolley
ST00497B	DNA strap belt w/integrated re-winding system
ST00498B	DNA strap thoracic belt w/re-winding system
ST00499B	STX499 – Thorax belts <small>Mandatory for 10G versions</small>
CB09028	Supporto fiss. tanker per Safe Bar

### Technical data <sup>1</sup>

Height <sup>(CM)</sup>	1430 mm
Width <sup>(CM/SM)</sup>	600 mm
Depth with all handles closed <sup>(CM)</sup>	1000 mm
Depth with all handles opened <sup>(CM)</sup>	1365 mm
Length with handles closed <sup>(SM)</sup>	1970 mm
Length with handles opened <sup>(SM)</sup>	2235 mm
Side bar dimensions	600x200 mm
Weight <sup>(With mattress)</sup>	29,5 ± 0,5 kg
Loading capacity	250 kg

<sup>1</sup>CM:Chair mode – SM: Stretcher mode

Class I Med.Device according to Dir. 93/42/CEE

CROSS CHAIR  
STRETCHER/CHAIR YELLOW/B  
CR00080C



CND Classification

V0899

Rev.2 (16/06/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)