

## CROSS CHAIR

### Trolley 10G for Stretcher/Chair



Trolley designed to be used with the Stretcher/Chair Cross Chair, for easier handling, loading and fixing inside ambulances. 10G certified if used with Spencer fastener R-Max.

#### Specific Characteristics

- Robust stainless steel frame
- Wheels  $\varnothing 200$  to increase sliding action
- Parking brakes on rear wheels
- Possibility to place in intermediate height to facilitate loading of Cross Chair
- Cross Chair fixing system easy to operate

#### Technical data

|   |                      |
|---|----------------------|
| Lenght  | 1980 $\pm$ 10 mm     |
| Width   | 570 $\pm$ 5 mm       |
| Sliding profile height (Raised trolley)                 | 850/900 mm           |
| Sliding profile height (Trolley in intermediate height) | 600/680 mm           |
| Sliding profile height (Lowered trolley)                | 280 mm               |
| Wheels  | $\varnothing$ 200 mm |
| Weight  | 35,5 $\pm$ 0,5 kg    |
| Loading capacity <sup>1</sup>                           | 150 kg               |

<sup>1</sup> Considering patient weight, equipment and accessories

#### Accessories

|          |   |
|----------|---|
| CR00080C | Cross Chair Stretcher/Chair                 |
| ST42702C | Fastener R-Max (Mandatory for 10g versions) |

Class I Med.Device according to Dir. 93/42/CEE

CROSS CHAIR TROLLEY 10G

CR00082C

CND Classification



V0899

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)