

## CROSS CHAIR

### 10G for Stretcher/Chair with trolley



Trolley with detachable Stretcher/Chair Cross Chair, for easier handling, loading and fixing inside ambulances. 10G certified if used with Spencer fastener R-Max.

#### Trolley Specific features

Robust stainless steel frame
Wheels $\varnothing$ 200 to increase sliding action
Parking brakes on rear wheels
Possibility to place in intermediate height to facilitate loading of Cross Chair
Cross Chair fixing system easy to operate

#### Chair stretcher Specific features

Configuration change takes less than 2 seconds
4-segment PVC Mattress with removable cover
Polyethylene board easy to clean
Backrest adjustable in 5 positions
Telescopic handles
Auto positioning rear handles
Suitable for use with Cross Chair Trolley
Light and robust
The side bars <i>Safe bar</i> contain the patient also in Chair configuration
Wheels $\varnothing$ 150 mm to increase sliding action
Trendelenburg antishock position

#### Accessories

IF01043B	IV holder TRACK 4
ST42702C	Fastener R-Max <sup>(Mandatory for 10g versions)</sup>

#### Trolley technical data

Lenght	1980 $\pm$ 10 mm
Width	570 $\pm$ 5 mm
Sliding profile height (Raised trolley)	850/900 mm
Sliding profile height (Trolley in intermediate height)	600/680 mm
Sliding profile height (Lowered trolley)	280 mm
Wheels	$\varnothing$ 200 mm
Weight	35,5 $\pm$ 0,5 kg
Loading capacity <sup>1</sup>	170 kg

<sup>1</sup> Considering patient weight, equipment and accessories

#### Chair stretcher technical data <sup>1</sup>

Height <sup>(CM)</sup>	1430 mm
Width <sup>(CM/SM)</sup>	600 mm
Depth with all handles closed <sup>(CM)</sup>	1000 mm
Depth with all handles opened <sup>(CM)</sup>	1365 mm
Length with handles closed <sup>(SM)</sup>	1970 mm
Length with handles opened <sup>(SM)</sup>	2235 mm
Side bar dimensions	600x200 mm
Weight <sup>(With mattress)</sup>	29,5 $\pm$ 0,5 kg

Class I Med.Device according to Regulation 2017/745/UE

CROSS CHAIR-10G SELF LOADING  
STRETCHER CHAIR  
CR00085

CND Classification



V0899

Rev.1 (28/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)