

CROSS TROLLEY

Self loading trolley for incubator



Self loading trolley with adjustable heights for transport of incubators. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier.

Specific Characteristics

The distance between incubator supports is adjustable in order to adapt to your own incubator

Integrated stainless steel instruments table is large and easy to clean

Color coded commands

2 fixed wheels and 2 pivoting ones of large diameter, increase sliding action

Brakes on rear wheels

Adjustable heights system

Foldable loading carriage allows movement in confined spaces easier

Accessories

ST42100A	Rear fastener for self loading stretchers
ST42200A	Front fastener for self loading stretchers

Technical data

Length	1980 ± 10 mm
Shortened length	1660 ± 10 mm
Width	580 ± 10 mm
Table dimensions	520x450 mm
Table height <small>(From main frame)</small>	165 ± 5 mm
Weight	34 ± 2 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height loading carriage wheel	From 52 to 70 cm
Frame material	Steel

Class I Med.Device according to Dir. 93/42/CEE

CROSS CART INCUBATOR HOLDER
F/CHROME

CR10256C

CND Classification V0899

Registration number 362328



Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it