

CROSS

Self Loading stretcher



Spencer Cross is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier.

Specific Characteristics

Polyethylene board easy to clean
 Color coded commands
 Adjustable backrest
 Antishock positions Trendelenburg
 2 fixed wheels and 2 pivoting ones of large diameter, increase sliding action
 Brakes on rear wheels
 Foldable side bars *Safe bar*
 Adjustable heights system
 Foldable loading carriage to save space

Accessories

ST42100A	Rear fastener for self loading stretchers
ST42200A	Front fastener for self loading stretchers
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
IF01047C	Track 4-30 telescopic IV holder
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

ST70000A	PVC anatomic mattress QMX/777
ST70002A	STX 702 two pieces black reflex metal belts

Technical data

Length	1980 ± 10 mm
Shortened length	1660 ± 10 mm
Width	580 ± 10 mm
Weight	45 ± 1 kg
Loading capacity	170 kg
Wheels	∅200 mm
Height loading carriage wheel	From 52 to 70 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

CROSS T P/BLACK F/YELLOW
 W/MAT. QMX777

CR50054E

CND Classification V0899

Registration number 88367



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it