

CROSS UP 8409

10G certified trolley



Self-loading trolley to be used with the stretcher Cross Up 8409. The legs fold away independently, activating the dedicated commands during loading procedures, and open automatically when unloading. The adjustable heights systems is easy to use, lifting and lowering procedures require minimal effort. 10G Certified if used with a certified Spencer fixing system.

Specific Characteristics

Color coded commands
 2 fixed wheels and 2 pivoting with large diameter to make movement easier
 Brakes on rear wheels
 Adjustable heights system to make easier the transfer procedures of the patient placed on the stretcher

Technical data

| | |
|------------------------------------|-----------------|
| Length | 1975 ± 10mm |
| Width | 580 ± 10mm |
| Stretcher support max. height | 910 ± 10mm |
| Stretcher support intermed. height | 680 ± 10mm |
| Stretcher support min. height | 280 ± 10mm |
| Height loading carriage wheel | 65 cm |
| Wheels | ø200 mm |
| Materials | Acciaio, Al, PE |
| Loading capacity | 250 kg |
| Weight | 31 ± 1 kg |

Accessories

| | |
|----------|---|
| ST42702C | R-Max fixing system <small>(Mandatory for 10G versions)</small> |
| CR00009E | STRETCHER S CROSS UP 8409 |
| CR00010E | STRETCHER T CROSS UP 8409 |
| CR00012E | STRETCHER TF CROSS UP 8409 |
| CR00013E | STRETCHER S CROSS UP 8409 PER B-BAK PIN |
| CR00014E | STRETCHER T CROSS UP 8409 PER B-BAK PIN |
| CR00015E | STRETCHER TF CROSS UP 8409 PER B-BAK PIN |

Class I Med.Device according to Dir. 93/42/CEE

CROSS UP TROLLEY 8409

CERTIFIED 10 G

CR60003E



CND Classification V0899

Registration number 792768

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it