

## CROSS UP 8409

10G certified trolley



Self-loading trolley to be used with the stretcher Cross Up 8409. The legs fold away independently, activating the dedicated commands during loading procedures, and open automatically when unloading. The adjustable heights systems is easy to use, lifting and lowering procedures require minimal effort. 10G Certified if used with a certified Spencer fixing system.

### Specific Characteristics

- Color coded commands
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Adjustable heights system to make easier the transfer procedures of the patient placed on the stretcher

### Technical data

Length	1975 ± 10mm
Width	580 ± 10mm
Stretcher support max. height	910 ± 10mm
Stretcher support intermed. height	680 ± 10mm
Stretcher support min. height	280 ± 10mm
Height loading carriage wheel	65 cm
Wheels	ø200 mm
Materials	Acciaio, Al, PE
Loading capacity	250 kg
Weight	31 ± 1 kg

### Accessories

ST42702C	R-Max fixing system <small>(Mandatory for 10G versions)</small>
CR00009E	STRETCHER S CROSS UP 8409
CR00010E	STRETCHER T CROSS UP 8409
CR00012E	STRETCHER TF CROSS UP 8409
CR00013E	STRETCHER S CROSS UP 8409 PER B-BAK PIN
CR00014E	STRETCHER T CROSS UP 8409 PER B-BAK PIN
CR00015E	STRETCHER TF CROSS UP 8409 PER B-BAK PIN

Class I Med.Device according to Dir. 93/42/CEE

CROSS UP TROLLEY 8409  
CERTIFIED 10 G

**CR60003E**



**CND Classification** V0899

**Registration number** 792768

Rev.2 (03/03/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)