

CROSS UP 8409 TWIST 10G

Trolley with adjustable heights



Self-loading trolley to be used with the stretcher Cross Up 8409. The legs fold away independently, activating the dedicated commands during loading procedures, and open automatically when unloading. The adjustable heights systems is easy to use, lifting and lowering procedures require minimal effort. 10G certified if used with the dedicated fastener R-MAX or S-MAX

Specific Characteristics

- Color coded commands
- Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting all 4 wheels (360 degrees).
- Brakes on rear wheels
- Adjustable heights system to make easier the transfer procedures of the patient placed on the stretcher
- Quick coupling mechanism to fix the 8409 stretcher

Technical data


Length	1975 ± 10mm
Width	580 ± 10mm
Stretcher support max. height ¹	910 ± 10mm
Stretcher support intermed. height	680 ± 10mm
Stretcher support min. height	310 ± 10mm
Height loading carriage wheel	65 cm
Wheels	ø200 mm
Materials	Steel, Al, PE
Loading capacity	250 kg
Weight	31 ± 1 kg

¹The indicated heights are a mean value of the stretcher support. The patient floor with the stretcher on the trolley is about 70mm higher.

Accessories

ST42702C	Fastener R-Max <small>(Mandatory for 10G versions)</small>
ST42707B	Fastener S-Max <small>(Mandatory for 10G versions)</small>
CR00009E	STRETCHER S CROSS UP 8409
CR00010E	STRETCHER T CROSS UP 8409
CR00012E	STRETCHER TF CROSS UP 8409
CR00013E	STRETCHER S CROSS UP 8409 FOR B-BAK PIN
CR00014E	STRETCHER T CROSS UP 8409 FOR B-BAK PIN
CR00015E	STRETCHER TF CROSS UP 8409 FOR B-BAK PIN

Class I Med.Device according to Dir. 93/42/CEE

TROLLEY 8409 TWIST CERTIFIED 10 G	
CR60004E	
CND Classification	V0899
Registration number	793227

Rev.1 (23/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it